
BD Accuri[™] C6 and BD Stemflow[™] Pluripotent Stem Cell Kits: Characterization of stem cell pluripotency

This experimental data demonstrates the characterization of human pluripotent stem cell preparations using the BD AccuriTM C6 personal flow cytometer and established combinations of cell surface and intracellular pluripotency markers.

Experiment	H9 human embryonic stem cells (hESCs) were dissociated and analyzed for expression of key pluripotency markers (SSEA-3, SSEA-1, Tra-1-81, Sox2, Nanog, and Oct3/4, SSEA-4). Data was collected and analyzed on a BD Accuri C6 personal flow cytometer. All data shown was gated based on the light scatter properties of H9 hESCs.
Material	Reagents Panel A: BD Stemflow™ Human Pluripotent Stem Cell Sorting and Analysis Kit (Cat. No. 560461), including SSEA-3, SSEA-1, and Tra-1-81 markers
	Panel B: BD Stemflow™ Human Pluripotent Stem Cell Transcription Factor Analysis Kit (Cat. No. 560589), including Sox2, Nanog, and Oct3/4 markers Panel C: BD Stemflow™ Human and Mouse Pluripotent Stem Cell Kit (Cat. No. 560477), including SSEA-1 and SSEA-4, Oct3/4 markers
	Additional Material BD Matrigel™ hESC-qualified matrix (Cat. No. 354277) BD™ Accutase™ Cell Detachment Solution (Cat. No. 561527)



Data



Discussion

The data demonstrates that rapid and multiparametric flow cytometry using the BD Accuri C6 personal flow cytometer and various BD Stemflow Pluripotent Stem Cell kits is effective for characterization of pluripotent stem cell preparations.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Class 1 Laser Product.

Accutase is a trademark of Innovative Cell Technologies, Inc. Alexa Fluor® is a registered trademark of Molecular Probes, Inc.

CyTM is a trademark of Amersham Biosciences Corp. Cy dyes are subject to proprietary rights of Amersham Biosciences Corp and Carnegie Mellon University and are made and sold under license from Amersham Biosciences Corp only for research and in vitro diagnostic use. Any other use requires a commercial sublicense from Amersham Biosciences Corp, 800 Centennial Avenue, Piscataway, NJ 08855-1327, USA.

