
BD Biosciences

BD FACSCanto™ System Consumables and Spares List



BD FACSCanto™ System Consumables and Spares List

DESCRIPTION	CAT. NO.
Air Filter	336303
Auxiliary Level Sensor 1-SW 24V 7.5 IN Stem	334911
Bal Seal (SIP)	343509
Barbed Fitting 1/16 IN inner diam. Tubing	343508
Coupling Insert Barb-elbow 1/4 IN Flow	333072
Filter 0.2 µm 440 cm2 33 mm Capsule PPL	336945
Filter 0.45 µm 50 mm D 3/16 1/4 Barbs PPL	336971
Filter Bypass Assembly	33576007
Fuse, 2.5A/250VAC, 5 x 20 mm	343565
O-Ring 0.1761 IN inner diam. x 0.318 IN outer diam.	340086
O-Ring 0.031 IN inner diam. x 0.44 IN outer diam.	340087
O-Ring, 0.926 IN inner diam. x 0.070 IN outer diam.	343616
Plunger Spring Eject	335512
Retainer Bal Seal 0.125 IN diam. SIT	336999
Retainer, Plunger	335516
Sheath Level Sensor 5-SW 24V 10.6 IN Stem, BD FACSCanto 6-Color	334914
Sheath Level Sensor 5-SW 120V 10.6 in Stem 6 Pin, BD FACSCanto system	338979
Spring Compression 0.480 IN outer x 0.035 IN diam.	335514
Waste Tank 10 Liter	340261
Tank, 5 Liter	333504
Test Tubes - 12 x 75 mm, 10/case (min. order 125)	343675
Tube, Silicone 1/32 IN inner x 1/8 IN outer diam. (50)	337033
Waste Level Sensor 6-SW 24V 10.6 IN Stem, BD FACSCanto 6-Color	334915
Waste Level Sensor 2-SW 120V 10.6in Stem 7 Pin, BD FACSCanto system	338978
Waste Tank Caps Pack of 12	33885407
Wet Cart Liquid Filter Assembly	33566307
Assembly Cover-Tube Loader	338709
Power Cord, 6-ft for USA (15A, 5-150/320-C13)	337219
Wire 0.005 IN diam. x 7.3 IN L 316LVM-SS	337592
HTS Kit Parts Replacement	644787
BD FACS Clean 5L	340345
BD FACSTFlow Sheath Fluid 20L	342003
BD FACS Shutdown Solution 5L	334224

Class 1 Laser Products.

For In Vitro Diagnostic Use (IVD).

BD, BD Logo, and all other trademarks are property of Becton, Dickinson and Company. © 2014 BD



BD Biosciences
bdbiosciences.com